

## CS162: Premium Matcha Shot - 8 sachets 6x(8x1g)

Uncontrolled if printed

		Version:	1H / 19.01.2023
		Label Code:	CS162-01-21
<b>General Information</b>			
Product Description:	Finely ground green tea leaf powder.		
Unit Size:	(8x1g)		
Country of Origin:	Japan		

**Organic Status**

Organic:	Yes
Cert. at Source:	JONA
Cert. at Point of Sale:	Soil Association (UK), SKAL (EU/ non-EU)

**Ingredients**

Organic Matcha green tea powder (100%).

**Allergen Information**

Cereals containing gluten*:	Absent
Peanuts*:	Absent
Nuts*:	Absent
Soybeans*:	Absent
Celery*:	Absent
Mustard*:	Absent
Sesame*:	Absent
Sulfur dioxide and sulphites*+:	Absent
Crustaceans*:	Absent
Eggs*:	Absent
Fish*:	Absent
Milk*:	Absent
Lupin*:	Absent
Molluscs*:	Absent
* and products thereof +(>10mg/kg or 10mg/l)	

**Nutritional Information**

Typical figures in g/100g	
Energy kJ	1330 kJ
Energy kCal	317 kcal
Fat	0.8 g
of which saturates	0.13 g
Carbohydrate	38.7 g
of which sugars	10 g
Protein	23.4 g
Salt	< 0.1 g

Nutritional values for the product as sold

#### Production Process

Production Process:	Harvest -> steam -> cool -> dry -> selecting -> Tencha Aracha  Receiving Tencha Aracha (unrefined tea) -> quality check -> grinding (magnet: 10000 gauss) -> selecting -> drying (at 80-100C. Magnet: 10000 gauss) -> magnet (10000 gauss) -> metal detection -> X ray scanner -> milling (magnet: 10000 gauss) -> sieving (60 mesh) -> packaging -> weighing -> X ray scanner (bulk only) -> outer packaging -> shipping
Pasteurisation Process:	No
Processing Aids:	No

#### Dietary Selling Points

Gluten Free:	No
Raw Food:	No
Kosher:	Yes
Free of Added Alcohol:	Yes
Yeast Free:	Yes
Non-GM:	Yes
Vegan:	Yes
Free from added sugar:	Yes
Free from palm oil:	Yes
Free from hydrogenated oils:	Yes
Free from artificial additives:	Yes

#### Additional Information

n/a

#### Storage & Shelf Life




Min. Shelf life from Production:	24	Does not exclude shelf-life lost in transport from manufacturer
Min. Shelf life on Delivery:	Please check with your Account Manager	
Recommended Storage:	Store in a cool, dry place.	
Delivery Conditions:	Ambient temperature (between 10°C and 25°C)	
After opening:	Once opened, reseal and refrigerate.	
Notes and Cautions:	n/a	
Usage Instructions:	Sachet: add 1 sachet to 30ml hot water and whisk well.  Matcha powder: add 1/4 teaspoon (1g) to 30ml of hot water or a drink of your choice and whisk well.	
Traceability Coding:		
Use By Date Location:		

#### Packaging Information

Trade Case Content:	6
Unit Size (g or ml):	(8x1g)
Unit Description:	8 sachets in paper box
Trade Case description:	Paper box
Shelf Ready Packaging:	Yes
Trade Case Configuration (Front x Side):	02 x 03

### CS162: Premium Matcha Shot - 8 sachets 6x(8x1g)

Uncontrolled if printed

 Unit		 Case		 Pallet	
Barcode:	5021554003236	Barcode:	05021554003229	<b>Block</b> (1000mmx1200mm) Layers/pallet 12 Cases/layer 70	
Net Weight (g):	8	Net weight (g):	48		
Gross Weight (g):	24.4000000953674	Gross weight (g):	179.400000572205		
Drained Weight (g):	0				
Depth (mm):	53	Depth (mm):	164	<b>EURO</b> (800mmx1200mm) Layers/pallet 12 Cases/layer 56	
Width (mm):	49	Width (mm):	100		
Height (mm):	124	Height (mm):	128		
Paper (g):	10	Paper (g):	33		
Glass (g):	0	Glass (g):	0		
Aluminium (g):	0	Aluminium (g):	0		
Steel (g):	0	Steel (g):	0		
Plastic (g):	6.4	Plastic (g):	0		
Type of Plastic:	OTHER	Type of Plastic:			
Wood (g):	0	Wood (g):	0		
Other (g):	0	Other (g):	0		
Type of Other:	PET/PE/MAL/PE	Type of Other:			

Languages on Label:	English	Sizes Available:	139D	4000 X 1g
	Spanish		CS159	10 X 40g
	Italian		CS160	1 X 1kg
	Portuguese		CS161	30 X 1g
	French			
	Dutch			
	Norwegian			
	Danish			
	Swedish			
	Finnish			

Marketing Information	
Suggestions for use:	Sachet: add 1 sachet to 30ml hot water and whisk well.  Matcha powder: add 1/4 teaspoon (1g) to 30ml of hot water or a drink of your choice and whisk well.

Pack Shot Image:

This document will be valid only on the date it is printed. For the most recent version please visit the Clearspring Dropbox website: <http://bit.ly/ClearspringUkLibrary> or <http://bit.ly/ClearspringExportLibrary>