

TEA0201.090A: Sencha Green Tea - loose 6x90g

Uncontrolled if printed

		Version:	2A / 14.07.2022
		Label Code:	TEA0201.090A.03
General Information			
Product Description:	Clearspring Organic Japanese Sencha is grown in the luscious rolling hills around Kyoto and Kyushu, areas known for their ideal climatic and soil conditions which produce some of the finest tea in Japan. The tender young leaves are steamed, rolled and dried immediately after harvesting to ensure freshness and aroma are locked in to bring you an authentic, flavourful Japanese Sencha tea experience.		
Unit Size:	90g		
Country of Origin:	Japan		

Organic Status

Organic:	Yes
Cert. at Source:	JONA
Cert. at Point of Sale:	Soil Association

Ingredients

Organically grown green tea leaves

Allergen Information

Cereals containing gluten*:	Absent
Peanuts*:	Absent
Nuts*:	Absent
Soybeans*:	Absent
Celery*:	Absent
Mustard*:	Absent
Sesame*:	Absent
Sulfur dioxide and sulphites*+:	Absent
Crustaceans*:	Absent
Eggs*:	Absent
Fish*:	Absent
Milk*:	Absent
Lupin*:	Absent
Molluscs*:	Absent
* and products thereof +(>10mg/kg or 10mg/l)	

Nutritional Information

Typical figures in g/100g	
---------------------------	--

Nutritional values for the product

Production Process

Production Process:	harvest -> steam -> cool -> rub and dry repeatedly -> dry -> select -> Aracha -> quality check -> selecting -> drying (microwave far infrared dryer upper heat 80-95°C, lower heat 90-130°C, 10-15min) -> blending (magnet: 10000 gauss) -> magnet: 12000 gauss -> metal detecting (Fe: 0.8mm/ AL: 1.5mm) -> X-ray scanner (SUS:1.5mm/ Ceramic: 3.0mm / Glass: 3.0mm/ Rubber: 4.0mm) -> Weighing -> Packaging
Pasteurisation Process:	N/A
Processing Aids:	N/A

Dietary Selling Points

Gluten Free:	No
Raw Food:	No
Kosher:	No
Free of Added Alcohol:	Yes
Yeast Free:	Yes
Non-GM:	Yes
Vegan:	Yes
Free from added sugar:	Yes
Free from palm oil:	Yes
Free from hydrogenated oils:	Yes
Free from artificial additives:	Yes

Additional Information




n/a

Storage & Shelf Life

Min. Shelf life from Production:	24	Does not exclude shelf-life lost in transport from manufacturer
Min. Shelf life on Delivery:	Please check with your Account Manager	
Recommended Storage:	Store in a cool, dark place.	
Delivery Conditions:	Ambient temperature (between 10°C and 25°C)	
After opening:	Reseal pack after use.	
Notes and Cautions:		
Usage Instructions:	Place 1 heaped teaspoon of tea per person in a teapot. Boil fresh water and leave to cool for a minute before pouring over the tea. Infuse for only 1-2 minutes and enjoy.	
Traceability Coding:		
Use By Date Location:		

Packaging Information

Trade Case Content:	6
Unit Size (g or ml):	90g
Unit Description:	Plastic Flexible Pouch with ziplock opening
Trade Case description:	Cardboard box with a label
Shelf Ready Packaging:	No
Trade Case Configuration (Front x Side):	01 x 06

 Unit		 Case		 Pallet	
Barcode:	5021554000952	Barcode:	05021554000938	Block (1000mmx1200mm) Layers/pallet 3 Cases/layer 24	
Net Weight (g):	90	Net weight (g):	540		
Gross Weight (g):	98	Gross weight (g):	769		
Drained Weight (g):	0				
Depth (mm):	82	Depth (mm):	245	EURO (800mmx1200mm) Layers/pallet 3 Cases/layer 24	
Width (mm):	150	Width (mm):	240		
Height (mm):	210	Height (mm):	120		
Paper (g):	0	Paper (g):	181		
Glass (g):	0	Glass (g):	0		
Aluminium (g):	0	Aluminium (g):	0		
Steel (g):	0	Steel (g):	0		
Plastic (g):	8	Plastic (g):	0		
Type of Plastic:		Type of Plastic:			
Wood (g):	0	Wood (g):	0		
Other (g):	0	Other (g):	0		
Type of Other:	(Matte-ink/ PET12/ ink/ adhesive/ ALPET12/ adhesive/ antistatic- LLDPE60µm	Type of Other:			

Languages on Label:		Sizes Available:	TEA020	6 X 90g
			1.090-GB	
			TEA080	1 X 20kg
			1.X20A	

Marketing Information

Suggestions for use:	
	Place 1 heaped teaspoon of tea per person in a teapot. Boil fresh water and leave to cool for a minute before pouring over the tea. Infuse for only 1-2 minutes and enjoy.

Pack Shot Image:



This document will be valid only on the date it is printed. For the most recent version please visit the Clearspring Dropbox website: <http://bit.ly/ClearspringUkLibrary> or <http://bit.ly/ClearspringExportLibrary>