3.6.3 Clearspring Customer Specification

TEA0906.085A: Matcha Sencha, Green Tea with Matcha - loose

Uncontrolled if printed

	Version	n: 2B / 22.05.2023
General Information	Label Code	E: TEA0906.085A.02
Product Description:	Clearspring Organic Matcha Sencha is a harmonious blend of Japanese Sencha, married with distinctly rich and exquisite ceremonial Matcha. Specially created by our Japanese expert tea producers using only the finest Organic green tea leaves, it delivers the perfect balance of flavours being both uplifting and soothing, and can be enjoyed at any time of the day.	
Unit Size:	85g	
Country of Origin:	Japan	

Organic Status	
Organic:	Yes
Cert. at Source:	JONA
Cert. at Point of Sale:	Soil Association

Ingredients

Sencha green tea leaves* (95%), matcha green tea powder* (5%) *organically grown

Allergen Information	
Cereals containing gluten*:	Absent
Peanuts*:	Absent
Nuts*:	Absent
Soybeans*:	Absent
Celery*:	Absent
Mustard*:	Absent
Sesame*:	Absent
Sulfur dioxide and sulphites*+:	Absent
Crustaceans*:	Absent
Eggs*:	Absent
Fish*:	Absent
Milk*:	Absent
Lupin*:	Absent
Molluscs*:	Absent
* and products thereof +(>10mg/kg or 10mg/l)	

Nutritional Information

N/A - product exempt under Reg. (EU) 1169/2011

3.6.3 Clearspring Customer Specification

TEA0906.085A: Matcha Sencha, Green Tea with Matcha - loose

Uncontrolled if printed

Production Proces	SS
Production Process:	SENCHA harvest -> steam -> cool -> rub and dry repeatedly -> dry -> select -> Aracha -> quality check -> selecting -> drying (microwave far infrared dryer upper heat 80-95°C, lower heat 90-130°C, 10-15min) MATCHA harvest -> steam -> cool -> rub and dry repeatedly -> dry -> select -> Aracha -> quality check -> crushing (magnet: 10000 gauss) -> selecting -> drying (microwave far infrared dryer upper heat 80-95°C, lower heat 90-130°C, 10-15min) -> magnet (10000 gauss) -> magnet (12000 gauss) -> metal detecting (fe:0.8mm/ AL: 1.5mm) -> X-ray scanner (SUS:1.5mm/ Ceramic: 3.0mm/ Glass: 3.0mm/ Rubber: 4.0mm) -> grinding (bead mill) (magnet: 10000 gauss) -> sieving (60 mesh) MIX SENCHA & MATCHA blending (magnet: 10000 gauss) -> magnet: 12000 gauss -> metal detecting (Fe: 0.8mm/ AL: 1.5mm) -> X-ray scanner (SUS:1.5mm/ Ceramic: 3.0mm/ Glass: 3.0mm/ Rubber: 4.0mm) -> Weighing ->
	Packaging
Pasteurisation Process:	n/a
Processing Aids:	n/a

Dietary Selling Points	
Gluten Free:	No
Raw Food:	No
Kosher:	No
Free of Added Alcohol:	Yes
Yeast Free:	Yes
Non-GM:	Yes
Vegan:	Yes
Free from added sugar:	Yes
Free from palm oil:	Yes
Free from hydrogenated oils:	Yes
Free from artificial additives:	Yes

Additional Information wa

Storage & Shelf Life			
Min. Shelf life from Production:	24	Does not exclude shelf-life lost in transport from manufacturer	
Min. Shelf life on Delivery:	Please chec	k with your Account Manager	
Recommended Storage:	Store in a cool, dry place.		
Delivery Conditions:	Ambient temperature (between 10°C and 25°C)		
After opening:	Store in a cool, dry place.		
Notes and Cautions:			
Usage Instructions:	Place 1 heaped teaspoon of tea per person in a teapot. Boil fresh water and leave to cool for a minute before pouring over the tea. For the best flavour infuse for 1-2 minutes and serve.		
Traceability Coding:			
Use By Date Location:			

Packaging Information		
Trade Case Content:	6	
Unit Size (g or ml):	85g	
Unit Description:	Round tin box with plug in lid with temper proof	
Trade Case description:	Cardboard with tin dividers	
Shelf Ready Packaging:	No	
Trade Case Configuration (Front x Side):	02 x 03	

3.6.3 Clearspring Customer Specification

TEA0906.085A: Matcha Sencha, Green Tea with Matcha - loose

Uncontrolled if printed

	Unit		Case	Pallet
Barcode:	5021554003984	Barcode:	05021554003991	Block
Net Weight (g):	85	Net weight (g):	510	DIOCK
Gross Weight (g):	158.10000000149	Gross weight (g):	1078.60000000894	(1000mmx1200mm)
Drained Weight (g):	0			Layers/pallet 3
Depth (mm):	64	Depth (mm):	212	Cases/layer 72
Width (mm):	64	Width (mm):	142	
Height (mm):	125	Height (mm):	129	TUDO.
Paper (g):	0	Paper (g):	130	EURO
Glass (g):	0	Glass (g):	0	(800mmx1200mm)
Aluminium (g):	0	Aluminium (g):	0	Layers/pallet 3
Steel (g):	73	Steel (g):	0	Cases/layer 36
Plastic (g):	0.1	Plastic (g):	0	
Type of Plastic:		Type of Plastic:		
Wood (g):	0	Wood (g):	0	
Other (g):	0	Other (g):	0	
Type of Other:		Type of Other:		

Languages on Label:	English Spanish	Sizes Available:	TEA080 6.X20A	1 X 20kg
	French		TEA090 6.085A-	6 X 85g
	Dutch Norwegian		GB	
	Danish			
	Swedish			

Marketing Information		
Suggestions for use:	Place 1 heaped teaspoon of tea per person in a teapot. Boil fresh water and leave to cool for a minute before pouring over the tea. For the best flavour infuse for 1-2 minutes and serve.	

Pack Shot Image:



This document will be valid only on the date it is printed. For the most recent version please visit the Clearspring Dropbox website: http://bit.ly/ClearspringUkLibrary or http://bit.ly/ClearspringExportLibrary