Clearspring

#### 3.6.3 Clearspring Customer Specification

TEA0702.036A: Hojicha Rstd Green Tea - bags 4x20 bags (36g) Uncontrolled if printed

		Version:	1D / 08.03.2024
General Information		Label Code:	TEA0702.036A.02
Product Description:	Organic fine cut Hojicha roasted green tea		
Unit Size:	20 bags (36g)		
Country of Origin:	United Kingdom		

Organic Status	
Organic:	Yes
Cert. at Source:	JONA
Cert. at Point of Sale:	Soil Association (UK), SKAL (EU/ non-EU)

#### Ingredients

Organically grown roasted green tea leaves

Allergen Information	
Cereals containing gluten*:	Absent
Peanuts*:	Absent
Nuts*:	Absent
Soybeans*:	Absent
Celery*:	Absent
Mustard*:	Absent
Sesame*:	Absent
Sulfur dioxide and sulphites*+:	Absent
Crustaceans*:	Absent
Eggs*:	Absent
Fish*:	Absent
Milk*:	Absent
Lupin*:	Absent
Molluscs*:	Absent
* and products thereof +(>10mg/kg or 10mg/l)	

#### Nutritional Information

#### Typical figures in g/100g

Nutritional values for the product

## TEA0702.036A: Hojicha Rstd Green Tea - bags 4x20 bags (36g)

Uncontrolled if printed

<b>Production Process</b>	
Production Process:	Harvest -> steam -> cool -> rub and dry repeatedly -> dry -> select -> Aracha -> aging (minimum 3 months) -> quality check -> selecting -> cut and sieve (10-30 mesh) -> drying (gas dryer: 120C/248F, 10min) -> blend -> magnet (10,000 gauss) -> store -> roast (200C/392F, 30 min) -> magnet (12,000 gauss) -> metal detecting -> X-ray scanner -> weigh -> packaging *Receipt of bulk material -> sieve -> pack-> metal detection -> packaging *Repack process for retail size only
Pasteurisation Process:	N/A
Processing Aids:	N/A

### **Dietary Selling Points**

Gluten Free:	No
Raw Food:	No
Kosher:	Yes
Free of Added Alcohol:	Yes
Yeast Free:	Yes
Non-GM:	Yes
Vegan:	Yes
Free from added sugar:	Yes
Free from palm oil:	Yes
Free from hydrogenated oils:	Yes
Free from artificial additives:	Yes

## Additional Information

n/a

#### Storage & Shelf Life

Storage & Shell Life			
Min. Shelf life from Production:	24 Does not exclude shelf-life lost in transport from manufacturer		
Min. Shelf life on Delivery:	Please check with your Account Manager		
Recommended Storage:	Store in cool, dry place		
Delivery Conditions:	Ambient temperature (between 10°C and 25°C)		
After opening:	Once the tea envelope has been opened, use tea bag to prepare tea immediately.		
Notes and Cautions:			
Usage Instructions:	Boil fresh wa	ater and pour over the teabag. Infuse for 2 minutes only for best flavour.	
Traceability Coding:			
Use By Date Location:			

#### Packaging Information

Trade Case Content:	4	
Unit Size (g or ml):	20 bags (36g)	
Unit Description:	20 tea bags with envelope in a paper carton box	
Trade Case description:	4 carton boxes in a cardboard box	
Shelf Ready Packaging:	Yes	
Trade Case Configuration (Front x Side):	01 x 04	

# Clearspring

#### 3.6.3 Clearspring Customer Specification

## TEA0702.036A: Hojicha Rstd Green Tea - bags 4x20 bags (36g)

Uncontrolled if printed

	Unit		Case	Pallet
Barcode:	5021554001690	Barcode:	05021554001683	Block
Net Weight (g):	36	Net weight (g):	144	BIUCK
Gross Weight (g):	79.0300003066659	Gross weight (g):	401.219999700785	(1000mmx1200mm)
Drained Weight (g):				Layers/pallet 18
Depth (mm):	68	Depth (mm):	281	Cases/layer 28
Width (mm):	78	Width (mm):	81	
Height (mm):	139	Height (mm):	144	EURO
Paper (g):	42.93	Paper (g):	85.1	EURO
Glass (g):	0	Glass (g):	0	(800mmx1200mm)
Aluminium (g):	0	Aluminium (g):	0	Layers/pallet 18
Steel (g):	0	Steel (g):	0	Cases/layer 20
Plastic (g):	0	Plastic (g):	0	
Type of Plastic:		Type of Plastic:		
Wood (g):	0	Wood (g):	0	
Other (g):	0.1	Other (g):	0	
Type of Other:	Cotton string (actual weight: 0.02g)	Type of Other:		

Languages on Label:	English Spanish Italian Portuguese French Dutch Norwegian Danish Swedish Finnish German	Sizes Available: TEA050 1 X 20kg 2.X20A
	Arabic	

Marketing Information		
Suggestions for use:	Boil fresh water and pour over the teabag. Infuse for 2 minutes only for best flavour.	

#### Pack Shot Image:



This document will be valid only on the date it is printed. For the most recent version please visit the Clearspring Dropbox website: http://bit.ly/ClearspringUkLibrary or http://bit.ly/ClearspringExportLibrary